



CGNE-62-1(1)

#24  
attach

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re the Application of	)	
Comai, et al.	)	Examiner: P. Moody
Serial No. 07/985,742	)	Art Unit: 1804
Filed: December 4, 1992	)	
For: FIGWORT PLANT PROMOTER	)	DECLARATION UNDER
<u>AND USES</u>	)	<u>37 CFR 1.131</u>

Honorable Commissioner of  
Patents and Trademarks  
Washington, DC 20231

Dear Sir:

We, Luca Comai, Margaret P. Sanger and Stephen Daniel  
Daubert do hereby declare as follows:

1. We are the inventors of the subject application.  
The work represented in the attached notebook pages was  
conducted in reference to the above-identified patent  
application in the United States at least prior to November  
13, 1988.

2. Figwort mosaic virus 34S promoter construct pFWP-  
101 is described in the subject patent application.

3. Photocopies of relevant pages from Margaret P.  
Sanger's experimental notebook are attached hereto as Exhibit  
A.

4. Electroporation of plant protoplast with pFWP-101  
was conducted as shown on notebook pages 67-69 in Exhibit A.

5. Electroporated protoplasts were analyzed and GUS expression confirmed as demonstrated on notebook pages 69-71 in Exhibit A.

DECLARATION

We declare that all statements made herein of our own knowledge are true and that all statements made on information and belief are believed to be true; and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.

Date: 8/20/93

By:

Luca Comai  
Luca Comai

Date: 8/23/93

By:

Margaret P. Sanger  
Margaret P. Sanger

Date: 8/24/93

By:

Stephen Daniel Daubert  
Stephen Daniel Daubert

enclosures: Exhibit A

I Next need to prepare enough DNA to test promoter region by transient expression system using GUS activity

inoculated 50 ml L.B (Con 5ug/ml) with pCGN-7304-101


- inoculated - 7304-101 (50ml) & extracted plasmid DNA

- Meanwhile to Sal I cut FMV-X2 & FMV-X4  
- phenol/chloroform extracted

- EtOH pptd.

- dried & resuspended in 20ul for next infection. Test on *Sativa stramonium*

- II
- ① FMV (Sal I cut) X2 - 1ul
  - ② FMV X4 (Sal I cut) - 1ul
  - ③ pCGN 7304-101 uncut - 1ul
  - ④ pCGN 860 (35S-HM-trn3') 0.5ul
  - ⑤ pCGN 852 (35S-LM-trn3') 0.5ul
  - ⑥ pCGN 844 (trn3'-HL-0.53') 0.5ul
  - ⑦ 7304 (A) - 1ul
  - ⑧ 7304 (B) - 1ul



Checked concentration  
via A260

CGN 7304-101	250 x dil = .273	---	2.7mg/ml
CGN 7304	250 x dil = .914	---	9.14mg/ml
CGN 3000	2500 x dil = .512	---	50mg/ml
CGN 860	1000 x dil = .337	---	13.5mg/ml
CGN 852	500 x dil = .274	---	54.6mg/ml

Handwritten signature/initials.

Need to check CGN 7000 in gel

- Electroporation of FMV promoter  $\beta$ -glucuronidase-GUS into tobacco protoplasts:

- Want to test out FMV full length RNA promoter region - to see if altered <sup>works</sup> ~~transcription~~  
 - If AUG's cause problems?

→ Generally - is there any promoter activity

Will test the following constructs for comparison w/ FMV:

- (A) CGN 7000 (mas 5' - GUS - mas 3') -- 50  $\mu$ g + 175  $\mu$ g pUC19 DNA  
 (→ may have combined into (A) by mistake - array will tell)
- (B) CGN 7304 (double 35S - GUS - mas 3' -- 50  $\mu$ g + 175  $\mu$ g pUC19
- (C) CGN 7304-101 (FMV - GUS - mas 3') -- 50  $\mu$ g + 175  $\mu$ g pUC19
- (D) pUC19 - no GUS 225  $\mu$ g DNA.

The protocol used was that of Gubow, Jones & is outlined in my lab book #3208-002-859+60

In brief:

- 20 Kanthi leaves / 30 ml Enzyme solution
- 17 Kanthi leaves / 30 ml Enzyme "
- infiltrated at 300 mbar ~ 10:40
- incubated in dark ~ 2 hrs → 1:15 P
- The material was agitated by running up & down wide bore 10ml Japanese pipettes & then narrower (normal) bore Japanese pipettes
- strained through 5  $\mu$ m filter
- protos were centrifuged out of Enzyme ~ 7 PM
- & washed 2X (#2 setting IEC clinical cent for 4.5 min each spin)
- counted protos - counted 3 sectors & averaged

protoplast count @ 416

(B) 415

(C) 529

$$\bar{x} = 453 (\times 10^4 / \text{ml})$$

∴ have 8mls at  $4.5 \times 10^6$  protos/ml

$$\approx 36 \times 10^6 \text{ protos}$$

If electroporate 9 samples at 3.5 million/ml  
then need  $9 \times 3.5 = 31.5 \times 10^6$  protos.

- collected, spun → suspended in ELECTROPORATION BUFFER
- Added DNAs & → 1ml/sample
- Electroporated ~ 60msc.  
~ 1200  $\mu$ Farads  
~ 225  $\mu$ g Total DNA/sample

Following electroporation the 1ml samples  
were added to 10ml culture medium in  
plant culture petri dishes (25x100mm)

- Put to incubate in dark @ 25°C for ~ 2 days.

checked protos in 2m - still looked  
OK.

- Proto's - O.K. - no microbial contamination.
- protos still intact & looking O.K.

So - need to harvest, extract, & assay for  
GUS activity.

collected protos by pipetting them into 15ml  
screw-cap centrifuge tube, running <sup>upside down</sup>  
3mls wash buffer → tube & centrifuging  
for 8min at  $\frac{1}{2}$  full (~500g) for 8min.  
(note had to respin 1 sample because disturbed  
pellet & found 4min at  $\frac{3}{4}$  full is also O.K.).

- ① - pipetted off supernatant & resuspended protos in 1ml Extraction Buffer;

### GUS Extraction Buffer \*

50mM  $\text{Na}_2\text{HPO}_4$  (pH 7.0) or  $\text{NaH}_2\text{PO}_4$

10mM BME (14.4mM  $\text{Na}_2\text{S}_2\text{O}_5$  = 0.7mM/l)

10mM  $\text{Na}_2\text{EDTA}$

0.1% Sodium Lauryl Sarcosine

0.1% Triton X-100

→ Plus 1mM PMSF  $\text{MW} = 174.2$  should add  $\sim 100\text{mM}$   $\sim 20\text{mg/l}$

- put in sonic bath for 7 x 60sec & checked for protoplast leakage -- poor -- tried some more still at least 30% intact  
= polytroned each sample for ~ 5 sec each  
- most protos were broken & rest looked damaged.

- ② - put into 1.5ml conical cent. tubes & n'uged for 4min in cold. Used supernatant as crude enzyme solution

### GUS Assay

- ← used 100ul enzyme prep  
← add 400ul of Assay buffer - Extraction buffer  
by MUG 1.25mM  
T 500ul  
100mM  
125ul 100mM MUG  
→ 10ml ext. buff.  
→ at 37°  
so final assay conc of 500ul  
vol = 1mM  
for 10ml/wt 125ul of a  
100mM MUG solution

incubate for various times

①

MUG = methylumbelliferone glucuronide

(4-methyl umbelliferyl  $\beta$ -D-glucuronide  $\text{MW} \sim 440$   $\sim 32.5$ )

\* Richard A Jefferson (1987) "Assaying chimeric genes in plants: the GUS gene fusion system." Plant molec. Biol. Reporter

# STOP BUTIR

- ① Stop reaction w/ 900ul of 0.2M  $\text{Na}_2\text{CO}_3$  (21.2g/l)  
 - this stops rx as well as fluorescence of mCherry umbelliferone

To correlate relative fluorescence with [MU]  
 set-up standards (excit = 365, emiss = 455)

Autobug recommends following Wavelength:

100 $\mu\text{M}$	= 0.1 $\mu\text{M}$	(5ul 10 $\mu\text{M}$ + 495 $\rightarrow$ 500N)	- 5 + 495
500 $\mu\text{M}$	0.5 $\mu\text{M}$	(25ul 10 $\mu\text{M}$ + 495 + 495)	25 + 495
1.0		(5ul 100 $\mu\text{M}$ + 495 EB)	50 + 450
2.0		(10ul 100 $\mu\text{M}$ + 490)	100 + 400
4.0		(20ul 100 $\mu\text{M}$ + 480)	200 + 300
8.0		(40ul 100 $\mu\text{M}$ + 460)	400 + 100


Added substrate to extracts at 3:30 PM

- took first time pts @ 1 hr 15 min
- took 2nd time pt. @ 2 hr 45 min

1st looked at samples on UV transilluminator

- (A) = CGN 7000 + CGN 7304 50+50 + 175  $\mu\text{L}$  DNA
- (B) = mixed pmt tip of CGN 7304 (is mainly no DNA) (PIS PLC)
- (C) = CGN 7304-10 (50ul + 175  $\mu\text{L}$  PCR)
- (D) = pUC19 (225  $\mu\text{g}$ )

Thus FMV-promoter region is active and may be 30 to a level similar to a single 35S of CaMV



Fluorometric assay  $\rightarrow$

